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On the right, a detailed report view is shown. It includes patient information (キビタンID, 患者ID, 患者氏名, 性別, 生年月日), examination details (依頼元施設, 検査日付, モダリティ, 検査部位), and a list of titles (タイトル) for the DICOM images. At the bottom of this view, a 'Viewer 起動' button is highlighted with a red box, with a blue callout box pointing to it that says 'Viewer 起動ボタンを押下して下さい。' (Please click the Viewer activation button).

Other visible elements include a search bar at the top left, a '検索' button, and a '新規登録' button in the top right of the report view. The interface is in Japanese and uses a light blue color scheme.

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The image shows two screenshots of a DICOM viewer interface. The left screenshot displays a split view with a sagittal CT scan on the left and an axial CT scan on the right. A red arrow points from a thumbnail of the axial scan in the bottom-left corner to the main axial view. A blue arrow points from the main axial view to the right screenshot. The right screenshot shows the same sagittal view on the left, but the right panel is now a larger, full-screen view of the selected axial slice, highlighted with a red border. A blue text box at the bottom explains the action.

検査日時	部位	Im	Se
2017/09/06 10:47:36	CR CHEST	1	1
2017/09/06 15:41:42	CR CHEST	1	1
2017/09/01 09:49:05	CT	321	4
2015/02/26 14:11:06	MR LIVER	71	3

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